



MERIDIAN^{alpha} cranial perforators



MERIDIAN^{alpha}



Introduction

Meridian^{alpha} – a new cranial perforator for burr hole trepanation from adeor medical. Available as sterile, single-use version in two most common sizes. We offer a manufacturing standard second to none, resulting in an unmatched cutting efficiency. The cutting head profile design and the advanced release mechanism minimize potential risk of dura damage. The centre tip on the inner blade acts as an anti-skid feature upon start of the perforation. The result is safe and convenient cranial access.

Einführung

Meridian^{alpha} – ein neues Perforator-Design für die Burr-hole-Trepanation von adeor medical. Die Meridian^{alpha} Perforatoren sind als steril verpacktes Einmalprodukt erhältlich. Mit einem bisher unerreichten Qualitätsstandard erzielen die Meridian^{alpha} Perforatoren Trepanationsresultate der Extraklasse. Das Meridian^{alpha} Schneidenprofil schützt die Dura effektiv vor unbeabsichtigter Perforierung.

Single-Use Perforators

Made in Germany

Unique Benefits

- › Safe release mechanism
- › Forms Dura protecting bone pad
- › Suitable for even and uneven bone
- › Unmatched sharpness
- › Universal Hudson fitting
- › Durable cutting efficiency
- › Start & stop when needed
- › Restarts after stopping
- › Safe non-skid start
- › Available in two most common used sizes
- › Adult and pediatric versions (single-use)

Produktmerkmale

- › Sicherer Auslösemechanismus
- › Knochenplättchen als Duraschutz
- › Einsetzbar bei ebenen und unebenen Schädelknochen
- › Unübertroffene Schärfe
- › Universeller Hudsonschaft
- › Startet auch nach Unterbrechung
- › Rutschsicherer Anlauf
- › Verfügbar in zwei Standardgrößen
- › Standard und pädiatrische Versionen lieferbar (single-use)



Single-Use Perforators

For use on skulls as thin as 1.0 or 3.0 mm



Single-Use Perforators

For use on skulls as thin as 1.0 or 3.0 mm



MERIDIAN^{alpha} SINGLE-USE CRANIAL PERFORATOR

PER 7-11-S 7.0 mm / 11.0 mm (inner / outer Ø)

- > for use on skulls as thin as 3.0 mm
- > sterile, single-use, made in Germany
- > alternative to AcraCut® 200-243



MERIDIAN^{alpha} SINGLE-USE CRANIAL PERFORATOR

PER 11-14-S 11.0 mm / 14.0 mm (inner / outer Ø)

- > for use on skulls as thin as 3.0 mm
- > sterile, single-use, made in Germany
- > alternative to AcraCut® 200-271 or 200-241



MERIDIAN^{alpha} SINGLE-USE CRANIAL PERFORATOR

PER 7-11-SP 7.0 mm / 11.0 mm (inner / outer Ø)

- > for use on skulls as thin as 1.0 mm
- > sterile, single-use, made in Germany
- > alternative to AcraCut® 200-253



MERIDIAN^{alpha} SINGLE-USE CRANIAL PERFORATOR

PER 11-14-SP 11.0 mm / 14.0 mm (inner / outer Ø)

- > for use on skulls as thin as 1.0 mm
- > sterile, single-use, made in Germany
- > alternative to AcraCut® 200-251

adeor is a registered trademark of adeor medical AG.

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing device are reserved and require our written approval.

Copyright 2022 adeor medical AG
01-032022-D/E
Printed in Germany
March 2022

Codman is a trademark of Integra LifeSciences Inc.
AcraCut is a trademark of ACRA-CUT Inc.
adeor is not affiliated with any of the above.

adeor ist eine eingetragene Marke der adeor medical AG.

Alle Rechte der Verbreitung, Nachdruck (auch auszugsweise) oder Einspeicherung und Rückgewinnung in Datenverarbeitungsanlagen aller Art sind vorbehalten und bedürfen unserer schriftlichen Genehmigung.

Copyright 2022 adeor medical AG
01-032022-D/E
Gedruckt in Deutschland
März 2022

Codman ist eine Marke der Integra LifeSciences Inc.
AcraCut ist eine Marke der ACRA-CUT Inc.
adeor ist mit keiner der genannten Firmen affiliert.

HEADQUARTERS

adeor medical AG
Martinshof 5
D-83626 Valley
Germany

Phone +49 (0) 8024 47717 0
Fax +49 (0) 8024 47717 99
Email office@adeor.com
Web www.adeor.com



Printed on recycled paper.
Gedruckt auf Recyclingpapier.

